Myopia control effect of Repeated Low-Level Red-Light Therapy in Chinese children


A Randomized, Double-Blind, Controlled Clinical Trial

Jing Dong, MD | Zhuoting Zhu, MD, PhD | Haifeng Xu, MD | Mingguang He, MD, PhD

Purpose:

To assess the efficacy and safety of RLRL therapy in controlling myopia progression compared to a sham device with only 10% of the original power.

Method:



Results:

- 1. A total of 111 children were included in the analysis (n = 56 in the RLRL group and n = 55 in the sham device control group).
- 2. In myopic children, RLRL therapy in 100% power significantly reduced myopia progression over 6 months compared with those treated with a sham device of 10% original power.
- 3. We found potential efficacy in reducing myopic progression and axial elongation with the sham device which used 10% of the original device's power, compared with the historic SVS control group.

Scan QR code to read the full article.





To find out more about the Repeated Low-Level Red-Light Therapy available via the Eyerising Myproclear, get in touch with your local Eyerising International team today.

Myopia control effect of Repeated Low-Level Red-Light Therapy in Chinese children

A Randomized, Double-Blind, Controlled Clinical Trial

Jing Dong, MD | Zhuoting Zhu, MD, PhD | Haifeng Xu, MD | Mingguang He, MD, PhD

Clinical trial outcome:

Table 2 Change in Cycloplegic Spherical Equivalent Refraction and Axial Length from Baseline to 6 Months

Outcome	Intervention Group (n = 56)	Sham Device Control Group (n = 55)	Mean Difference (95% Confidence Interval)	<i>P</i> Value∗
Primary				
Change of SER, D	0.06 ± 0.30	-0.11 ± 0.33	0.17 (0.05–0.29)	0.003
Secondary				
Change of AL, mm	0.02 ± 0.11	0.13 ± 0.10	-0.11 (-0.15 to -0.07)	< 0.001

AL = axial length; D = diopter; SER = spherical equivalent refraction.

Data are presented as mean ± standard deviation, unless otherwise indicated. Bold indicates statistically significant.

* One sided.

